Levenhuk DTX 720 WiFi Microscope connection to iOS/Android device (Android smart phone taken as an example)


Fig. 1: Connect the device

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

Fig. 4: Enter the user name and password and get the camera listed

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EN

Levenhuk DTX 720 WiFi Microscope connection to iOS/Android device (sample: Android smart phone)

1. Power on the microscope by moving the slider to WiFi position;
2. Turn on the Wi-Fi of your iOS/Android device, search the Wi-Fi digital microscope from the Wi-Fi Device list. The Wi-Fi Microscope name is WIFICAMxxxx (where xxxx stands for four-digit random number), as follows (see Fig. 1)
3. The password to connect the device is **123456789**.
4. Start the installed Levenhuk DTX application, as follows (see Fig. 2):
5. Click **Add new camera**, then click **Search**.
The searched Wi-Fi camera will be shown on the Search result list. Click the camera name (for example, WIFICAM). You will also see the camera ID (for example, OBJ-002864-STBZD) (see Fig. 3).
6. Enter the user name as **admin**, the default password as **123456** (see Fig. 4). Then the online Wi-Fi camera is listed, and you can start working with the microscope.

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